



Rapid 6 FOAM

All-in-one cleaning & disinfection solution

Rapid 6 Foam is a powerful, water-based solution for the rapid cleaning and disinfection of hard surfaces. Specifically designed to work with existing cleaning practices to provide high grade sporicidal disinfection throughout medical settings.



Overview

Using a unique synergistic blend of surfactants, biocides, quaternary ammonium compounds and PHMB, Rapid 6 Foam provides targeted performance against:

- Gram-positive bacteria
- Gram-negative bacteria
- Fungi
- Mycobacteria
- Enveloped viruses
- Bacterial endospores

Independently tested using 'real-world' conditions and contact times, Rapid 6 Foam has also been tested and proven to provide 4 hour residual protection against bacteria & spores.

Rapid 6 is available in both a Foam and Wipe format - which can be used individually or in combination together. This reduces the number of products hospitals require and allows flexibility for hospital-specific cleaning practices.

Rapid 6 Foam provides a focused foaming spray which helps users clearly see the surface area where the foam is being applied. This helps ensure even and consistent surface coverage – which is important, as inadequate biocide application can lead to micro-organism cross resistance developing.

Using Rapid 6 Foam in combination with Rapid 6 Wipes, the user is able to:

- Extend the cleaning power and increase the performance of the wipe
- Reduce the possibility of bio-burden occurring
- Ensure uniform biocide application

Rapid 6 Foam has also been compatibility tested on a range of other wipe materials to demonstrate consistent and effective biocidal performance. These include detergent wipes and dry microfibre cloths.

Test results

Spores

Pathogen	Test	Pass time
<i>Clostridium difficile</i>	EN 13704	1 min.
<i>Bacillus cereus</i>	EN 14347	1 min.
<i>Bacillus subtilis</i>	EN 14347	1 min.

Fungi

<i>Aspergillus niger</i>	EN 1275	5 min.
	EN 1650	1 min.
	EN 13697	15 min.
	EN 13624	5 min.
<i>Candida albicans</i>	EN 14562	15 min.
	EN 1275	1 min.
	EN 1650	5 min.
	EN 13697	1 min.
	EN 13624	1 min.
	EN 14562	15 min.

Mycobacteria

<i>Mycobacterium avium</i>	EN 14348	5 min.
	EN 14563	5 min.
<i>Mycobacterium terrae</i>	EN 14348	5 min.
	EN 14563	5 min.

Bacteria

<i>Enterococcus hirae</i>	EN 1276	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	5 min.
<i>Escherichia coli</i>	EN 1276	1 min.
	EN 13697	5 min.
<i>Klebsiella pneumoniae</i>	EN 1276	1 min.
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	EN 1276	1 min.
	EN 1040	1 min.
	EN 1276	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
<i>Staphylococcus aureus</i>	EN 14561	5 min.
	EN 1040	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	5 min.

Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	EN 1276	1 min.
	EN 1276	1 min.
	EN 1276	1 min.

Viruses

H1N1 Swine Flu†	EN 14476†	1 min.
-----------------	-----------	--------

Other results

4 hour residuality testing against bacteria and spores
Practical "in use" testing against *C. difficile* spores

Contact us for further details on our test results. †Screen testing

RAPID 6 FOAM IS DESIGNED FOR PROFESSIONAL USE ONLY. USE DISINFECTANTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE. ALWAYS FOLLOW CLEANING PROCEDURES AND GUIDELINES.

* For further details, please contact us.

How to use

Rapid 6 Foam + a wipe of your choice



Ready to use - simply spray and wipe.



Use with Rapid 6 Wipes to boost performance by up to 100 times!*

OR



Also proven to provide effective cleaning and disinfection with microfibre cloths or Endurocide® Detergent Wipes.



Where to use

Rapid 6 Foam is designed for use on all hard surfaces (excluding medical devices) and can be used to clean & disinfect amongst others:

- Isolation rooms
- Laboratories
- Mortuaries
- Sluice rooms
- Blood spills
- Bodily fluid spills
- Macerators
- Beds & mattresses
- Dressing trolleys
- Walls & floors
- Steel gurneys
- Door handles
- Commodes
- Work stations
- Toilets, sinks & taps

Compatibility

Rapid 6 Foam is compatible with:

- Steel
- Cotton
- PVC
- Polyester
- Rubber
- Wood

Pack sizes

Description: trigger sprays	Case Size	Order Code	Description: refills	Case Size	Order Code
500ml foaming trigger spray	12 per case	R6-F-500	5 Litres refill	2 per case	R6-F-2X5
1 Litre foaming trigger spray	12 per case	R6-F-1000	5 Litres refill	4 per case	R6-F-4X5