

Rapid 6 Foam • Test Results

Wipes & Foam testing - *C.difficile* spores

Testing against *C.difficile* spores when Rapid 6 Foam is used with different wipe applications.

Method: 4 tiles were treated with *C.difficile* spores, then left to dry. The treated tiles were then cleaned using four different methods: Rapid 6 Wipes only; Rapid 6 Foam & microfibre cloths; Rapid 6 Foam & Endurocide Detergent Wipes; Rapid 6 Foam & Rapid 6 Wipes. Swabs were taken from both the cleaned tiles and the wipe after 1, 2 and 5 minutes.

Test conditions: done under dirty conditions.

Method used	Log. reduction	
	Surface	Wipe
Rapid 6 Foam + Rapid 6 Wipes		
1 minute	> 5.46	> 5.46
2 minutes	5.74	5.06
5 minutes	> 6.97*	> 6.97*
Rapid 6 Wipes only		
1 minute	> 3.47	> 3.24
2 minutes	4.87	4.46
5 minutes	> 6.41*	> 6.41*
Rapid 6 Foam + microfibre cloths		
1 minute	3.86	> 4.01
2 minutes	4.65	4.97
5 minutes	> 6.97*	> 6.27
Rapid 6 Foam + Endurocide Detergent Wipes		
1 minute	3.64	3.10
2 minutes	4.45	3.92
5 minutes	> 6.68	4.68

EN 17126 Spores

A quantitative suspension test designed to evaluate sporicidal activity in the medical area.

Pass: ≥ 4 log reduction in 15 minutes at 20°C under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Bacillus cereus</i>	✓	> 4.00	4 min.
<i>Bacillus subtilis</i>	✓	> 4.00	4 min.
<i>Bacillus subtilis</i>	✓	> 4.00	4 min.

EN 13704 Spores

A suspension test designed to evaluate sporicidal activity.

Pass: ≥ 3 log reduction in 60 minutes at 20°C under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Clostridium difficile</i>	✓	3.21	1 min.
	✓	4.06	2 min.
	✓	> 6.71*	5 min.

EN 13704 (liquid expelled from wipe) Spores

A suspension test designed to evaluate sporicidal activity.

Pass: ≥ 3 log reduction in 60 minutes at 20°C under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Clostridium difficile</i>	✓	3.86	1 min.
	✓	4.46	2 min.

EN 14347 Spores

A quantitative suspension test designed to evaluate basic sporicidal activity.

Pass: ≥ 3 log reduction in 30 minutes at 20°C.

Pathogen	Pass	Log. reduction	Test time
<i>Bacillus cereus</i>	✓	4.49	1 min.
<i>Bacillus subtilis</i>	✓	> 6.61*	1 min.

EN 14348 Mycobacteria

A suspension test designed to evaluate mycobactericidal activity in the medical area.

Pass: ≥ 4 log reduction in 5 minutes at 20°C under clean conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Mycobacterium avium</i>	✓	4.96	5 min.
<i>Mycobacterium terrae</i>	✓	5.01	5 min.

EN 14563 Mycobacteria

A quantitative carrier test designed to evaluate mycobactericidal activity of disinfectants in the medical area in the conditions in which they are used. **Pass:** ≥ 4 log reduction in 60 minutes under clean conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Mycobacterium avium</i>	✓	4.49	5 min.
<i>Mycobacterium terrae</i>	✓	4.11	5 min.

Other results

Spores - residuality tested	Tested against <i>C.difficile</i> spores - provides sporicidal activity 4 hours after application (log 3.8).
Bacteria - residuality tested	Tested against <i>E.coli</i> - provides antimicrobial activity 4 hours after application (log 5.55).
Bacteria & fungi - residuality tested	Tested against MRSA, <i>Salmonella</i> and <i>Candida albicans</i> - provides > 99.9% antimicrobial activity 3 hours after application

EN 1040 Bacteria

A suspension test designed to evaluate basic bactericidal activity.

Pass: ≥ 5 log reduction in 5 minutes with no interfering substance.

Pathogen	Pass	Log. reduction	Test time
<i>Pseudomonas aeruginosa</i>	✓	> 6.40*	1 min.
<i>Staphylococcus aureus</i>	✓	> 6.37*	1 min.

EN 1276 Bacteria

A suspension test designed to evaluate bactericidal activity.

Pass: ≥ 5 log reduction in 5 minutes under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Escherichia coli</i>	✓	> 6.28*	1 min.
<i>Enterococcus hirae</i>	✓	> 6.45*	1 min.
ESBL (Extended Spectrum Beta-Lactamase) <i>E.coli</i>	✓	> 6.41*	1 min.
<i>Klebsiella pneumoniae</i>	✓	> 6.43*	1 min.
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	✓	> 6.32*	1 min.
<i>Pseudomonas aeruginosa</i>	✓	> 6.32*	1 min.
<i>Salmonella enteritidis</i>	✓	> 6.46*	1 min.
<i>Salmonella typhimurium</i>	✓	> 6.32*	1 min.
Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	✓	> 6.48*	1 min.

EN 14561 Bacteria

A quantitative carrier test designed to evaluate bactericidal activity of disinfectants in the medical area in the conditions in which they are used. **Pass:** ≥ 5 log reduction in 60 minutes under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Enterococcus hirae</i>	✓	> 6.18*	5 min.
<i>Escherichia coli</i>	✓	> 6.48*	5 min.
<i>Pseudomonas aeruginosa</i>	✓	> 6.30*	5 min.
<i>Staphylococcus aureus</i>	✓	> 6.43*	5 min.

EN 13697 Bacteria

A quantitative surface test designed to evaluate bactericidal activity of disinfectants on non-porous surfaces in the medical area. **Pass:** ≥ 4 log reduction in 5 minutes under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Enterococcus hirae</i>	✓	6.00	5 min.
<i>Escherichia coli</i>	✓	> 6.23*	5 min.
<i>Pseudomonas aeruginosa</i>	✓	4.38	5 min.
<i>Staphylococcus aureus</i>	✓	> 6.46*	5 min.

EN 13727 Bacteria

A quantitative surface test designed to evaluate bactericidal activity of disinfectants on instruments in the medical area. **Pass:** ≥ 5 log reduction in 60 minutes under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Enterococcus hirae</i>	✓	> 6.26*	1 min.
<i>Pseudomonas aeruginosa</i>	✓	> 6.34*	1 min.
<i>Staphylococcus aureus</i>	✓	> 6.45*	1 min.

EN 1275 Fungi

A quantitative suspension test designed to evaluate basic fungicidal activity.

Pass: ≥ 4 log reduction in 15 minutes at 20°C with no interfering substance used.

Pathogen	Pass	Log. reduction	Test time
<i>Aspergillus niger</i>	✓	4.83	5 min.
<i>Candida albicans</i>	✓	> 5.51*	1 min.

EN 1650 Fungi

A suspension test designed to evaluate fungicidal activity.

Pass: ≥ 4 log reduction in 15 minutes at 20°C under clean conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Aspergillus niger</i>	✓	4.59	1 min.
<i>Candida albicans</i>	✓	5.11	5 min.

EN 13697 Fungi

A quantitative suspension test designed to evaluate fungicidal activity of disinfectants.

Pass: ≥ 3 log reduction in 15 minutes at 20°C under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Aspergillus niger</i>	✓	> 4.88	15 min.
<i>Candida albicans</i>	✓	> 5.61	1 min.

EN 13624 Fungi

A suspension test to evaluate fungicidal activity for instruments used in the medical area.

Pass: ≥ 4 log reduction in 60 minutes at 20°C under dirty conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Aspergillus niger</i>	✓	5.99	5 min.
<i>Candida albicans</i>	✓	> 5.46*	1 min.

EN 14562 Fungi

A quantitative carrier test designed to evaluate fungicidal activity in the medical area.

Pass: ≥ 4 log reduction in 60 minutes at 20°C under clean conditions.

Pathogen	Pass	Log. reduction	Test time
<i>Aspergillus niger</i>	✓	4.87	15 min.
<i>Candida albicans</i>	✓	4.25	15 min.

EN 14476 Viruses

A quantitative suspension test designed to evaluate virucidal activity of disinfectants in the medical environment.

Pass: ≥ 4 log reduction in 60 minutes at 20°C under clean conditions.

Pathogen	Pass	Log. reduction	Test time
H1N1 Influenza A Virus (Swine Flu)	✓	≥ 4.13	1 min.

*Minimum log. reduction - limited by culture